


Applications of DNA sequencing

Overview

- DNA Sequencing Applications
 1. Resequencing
 2. Capture/Targeted Resequencing
 3. RAD/ddRAD/GBS
 4. Bisulfite Sequencing
 5. De novo assembly
- Reference Genome availability
- N50

Whole Genome Shotgun Sequencing

- Start with genomic DNA
- DNA is sheared into fragments
 - Physical
 - Acoustic shearing (Covaris)
 - Sonication (Bioruptor)
 - Hydrodynamic force (Hydroshear)
 - Enzymatic (transposase, DNase I)
 - Chemical
- Ideally, would like a very uniform size selection
 - paired end: depends on kit, from 200-600bp
 - mate pairs: 3-20 Kbp



App 1: Whole Genome Resequencing

- Sequencing multiple individuals from the same species
- Reference genome is already available
- Discover variations in the genomes between and within samples
 - mutations
 - insertions
 - deletions
 - rearrangements
 - copy number changes

How long do the reads need to be?

For the human genome, estimates are:

25mers = 80% unique coverage

43mers = 90% unique coverage

But longer is better for repeats, more complex genomes.

Quality Assessment



Trimming



Quality Assessment



Mapping to a Reference



Visualization



Calling variants



Assessing functional impact of variants




Submit to SRA

Workflow

App 2: Sequence Capture/Targeted Resequencing


- Reduce costs by only sequencing regions of interest
 - Increase multiplexing while maintaining deep coverage
 - Different target enrichment strategies
 - PCR amplification of many regions (QIAGEN)
 - Array-based Hybridization (NimbleGen)
 - In solution Hybridization* (Agilent and NimbleGen)
- Most common use is exome capture.
 - Can also use to capture exomes plus UTRs and miRNAs
 - Can also target large continuous regions, such as QTLs (quantitative trait loci)



Exome sequencing with Human Nextera Rapid Capture Exomes Kit


- >340,000 95mer probes
- Genomic footprint of 62Mb (~2% of the genome)

Start with a normal DNA library.



Exome sequencing with Human Nextera Rapid Capture Exomes Kit

- >340,000 95mer probes
- Genomic footprint of 62Mb (~2% of the genome)



D. Enrich using streptavidin beads

Exome sequencing with Human Nextera Rapid Capture Exomes Kit

- >340,000 95mer probes
- Genomic footprint of 62Mb (~2% of the genome)




Problems

- Off target fragments get sequenced
- Coverage is not uniform
- Genetic heterogeneity – do the probes work in all individuals?

Current human kits:

- 80% of bases are covered at 10X with 4Gb of sequence



Performance comparison of four exome capture systems for deep sequencing. Chilamakuri et al 2014.

“For all the technologies, 25 million reads were sufficient to cover about 80% of target bases with at least 10× depth, with the exception of the Nextera technology, which covered only about 60% of target bases with the same number of reads”

Workflow

- Essentially the same as a whole genome resequence project
- Depending on variants of interest, may need to assess where coverage failed

App 3: GBS/RADSeq

Polymorphic marker- An essential genomic tool for:

- Population Structure
 - Association mapping
 - Pedigree mapping
 - QTL mapping
 - Phylogeny
- Use the high volume and low cost of sequencing to replace SNP chips and microsatellites
 - How to efficiently use NGS for discovering markers?
 - How to efficiently use NGS to do the genotyping?

Restriction Site-associated DNA sequencing (RADSeq)


- Developed Baird et al 2008
- Identify and score thousands of genetic markers
- Randomly distributed across the target genome
- From many individuals using Illumina technology
- Subsampling only at specific sites defined by restriction enzyme


Patented by University of Oregon, licensed by:



Non profit use at universities can be licensed for free

Baird NA, Etter PD, Atwood TS, Currey MC, Shiver AL, et al. (2008) Rapid SNP Discovery and Genetic Mapping Using Sequenced RAD Markers. PLoS ONE 3(10): e3376.






1° reads	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTT	TG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTT	TG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTT	TG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTT	TG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTT	TG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTT	TG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTT	TG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTT	TG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTT	TG			
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTTT	TG			
2° reads	TGC	G	GGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTT	TG	
	TGCAGGACACACAGGAGCTGAGCCATTTCCTGCGGGCTCC	A	GACCAAACGTT	C	TG		
	TGCAGGACACACAGGAGC	A	GAGCCATTTCCTGCGGGCTCC	C	GACCAAACGTTT	TG	
	TGCAGGACACACAGGAGCTGAGCCATTTC	T	TGCGGGCTCC	A	GACCAAACGTTT	TG	
	TGCAGGACA	G	ACAGGAGCTGAGCCATTTCCTGCGGGCTCC	C	GACCA	G	ACGTTT

ddRAD

- May have the problem of too many sites across the genome (even if you use a rare cutter)
- Need a way to more accurately control the number of loci sequenced
- Double digest RAD
- Peterson et al 2012
- Also patented




Peterson et al (2012) Double Digest RADseq: An Inexpensive Method for De Novo SNP Discovery and Genotyping in Model and Non-Model Species. PLoS ONE 7(5): e37135.

ddRAD

- Double digest RAD
- (Peterson et al 2012)
- Simpler and cheaper library construction
- restriction digest with two enzymes simultaneously
- eliminate random shearing and end repair
- explicitly use size selection
- Sequence fragments generated by cuts with both REs and which fall within the size-selection window


A

RAD sequencing



B

double digest RADseq




Peterson et al (2012) Double Digest RADseq: An Inexpensive Method for De Novo SNP Discovery and Genotyping in Model and Non-Model Species. PLoS ONE 7(5): e37135.

GBS: Genotyping by Sequencing

- Elshire et al 2011
- Increased efficiency and cost benefits
- Reduced sample handling
- Methylation-sensitive REs used to filter out the repetitive portion of the genome
- Better barcoding system
- Fewer steps
- Free to use and to sell (ie not licensed)

*Last year paid \$46 per sample
96-plex, 1 plate
Including SNP calling
Including enzyme optimization
External rate at Cornell



Workflow


- Informatics can be difficult if you don't have a genome
- Much easier if you have a reference genome – can be similar to whole genome resequencing
- If not, some software packages:
 - Universal Network Enabled Analysis Kit (UNEAK) – part of the TASSEL software suite
 - Stacks

 Stacks



App 4: Bisulfite Sequencing

- DNA methylation
 - First discovered epigenetic mark
 - methyl groups are added to DNA
 - Suppresses transcription
 - Adenine and Cytosine can be methylated in prokaryotes
 - Only cytosine is methylated in eukaryotes



Bisulfite treatment

- How to figure out where methylation occurs while sequencing?
- Treatment of DNA with bisulphite:
 - Unmethylated cytosine -> uracil
 - 5-methylcytosine stays the same
- Sequencing can yield single- nucleotide resolution of methylation patterns

Problems:

- Incomplete conversion
- DNA degradation during conversion
- 5-methylcytosine and 5-hydroxymethylcytosine both read as a C in bisulphite sequencing



Bisulfite conversion



PCR



Workflow

- Need special software for mapping
 - Bismark
 - BSMap
 - BSMapper
- Downstream analysis
 - MethyKit – statistics, visualization, tiling windows



Reference Genomes

All the methods we talked about so far depend on (or are easier with) a reference genome.

How many genomes are out there to use as a base for mapping reads?

UCSC Genome Browser

- Mammal (49)
- Other Vertebrate (24)
- Deuterostome (3)
- Insect (13)
- Nematode (6)
- Other (2) – yeast and sea hare
- Viruses (1) - Ebola


UCSC Genome Browser on Ebola virus Sierra Leone 2014 (G3683/KM034562.1/eboVir3) Assembly


move <<< << < > >> >>> zoom in 1.5x 3x 10x base zoom out 1.5x 3x 10x 100x

KM034562v1:1-18,957 18,957 bp.

[hg38 replaces hg19 as default human assembly](#)

KM034562v1 KM034562v1 NP VP35 VP40 KM034562v1 GP VP30 VP24 KM034562v1 L





Reference Genomes

NCBI Genome

NCBI Genome has 4 levels:

- Complete - all chromosomes are gapless and have no runs of 10 or more ambiguous bases (Ns), there are no unplaced or unlocalized scaffolds, and all the expected chromosomes are present
- Chromosome - there is sequence for one or more chromosomes, gaps OK.
- Scaffold - some sequence contigs have been connected across gaps to create scaffolds, but the scaffolds are all unplaced or unlocalized
- Contig - nothing is assembled beyond the level of sequence contigs

NCBI Genome Records

- Viruses
 - 5,673 genomes
 - 5,639 complete (> 99%)
- Prokaryotes
 - 73,708 genomes
 - 5894 complete (8.0%)
 - 1024 chromosome level (1.4%)
- Eukaryotes
 - 3,494 genomes
 - 22 complete (< 1%)
 - 440 chromosome level (12.6%)

Sequencing for *de novo* Assembly

- Reconstructing the original full DNA molecules from (short) read fragments
- Jigsaw puzzle
- How do the pieces fit together? (overlap)
- Missing pieces (sequencing bias)
- Dirty pieces (sequencing error, real biological variation)



An example

A small “genome”:

Friends,

Romans,

countrymen

lend me your ears;

Reads:

ds, Romans, count
ns, countrymen, le
Friends, Rom
send me your ears;
cryman, lend me

Overlaps:

Friends, Rom
ds, Romans, count
ns, countrymen, le
crymen, lend me
send me your ears

Consensus:

Friends, Romans, countrymen, lend me your ears;

Sequencing for *de novo* Assembly

- Strategy differs significantly from resequencing
- Spectrum of difficulty:
 - Size
 - Repetitiveness
 - Polyploidy
 - Heterozygosity


From 1990s-2005, genomes were sequenced with Sanger technology (7-10X coverage) and BAC-based physical maps

Why don't we just do it the way we used to do it?

- Too much money!

NGS Strategies

- New Strategy 1: Ultra high throughput Illumina + Variety of sequence libraries
 - Mate pairs at a range of distances: 5kb, 10kb, 20kb, 40kb
- New Strategy 2: PacBio + Illumina
 - Error correct the long PacBio reads with Illumina
- New Strategy 3: PacBio only
 - More expensive.
- Challenges:
 - Initial contig build is computationally intensive
 - Many assembly algorithms require 100s of Gb of RAM to Tbs of RAM



Strategy 2.

Workflow

1. Clean all contaminants out of reads (much more critical than for resequencing!)
2. Build contigs
(There are many strategies for building contigs – next lesson we will learn about those and about why contigs don't span the whole chromosome)

Reads	ATGGCATTGCAA
	TGGCATTGCAATTG
	AGATGGTATTG
	GATGGCATTGCAA
	GCATTGCAATTGAC
	ATGGCATTGCAATT
	AGATGGTATTGCAATTG
Consensus Sequence	AGATGGCATTGCAATTGAC


<http://gcat.davidson.edu/phast/>


Workflow (cont)

3. Scaffold

Mate pairs

....20kb....





Mate pairs order and orient the contigs against each other.
Still leave gaps.

Workflow (cont)

4. Gap fill

- Use your existing reads
- Sanger sequencing

5. Assess quality and internal consistency

- Expected size?
- How many reads map consistently?

N50 statistic

- the length for which the collection of all contigs of that length or longer contains at least 50%
- similar to a mean or median of contig lengths
- used widely in genome assembly, especially in reference to contig lengths within a draft assembly.

N50 Example

- Assembly A contains six contigs of lengths:
 - 80 kbp, 70 kbp, 50 kbp, 40 kbp, 30 kbp, and 20 kbp
 - Sum size of assembly A is 290 kbp
 - N50 contig length is 70 kbp
 - “Half of the assembly is contained in contigs of 70kbp or greater”
 - If you randomly selected a location, 50% of the time it would be in a contig of 70kbp or greater
- Assembly B contains eight contigs of lengths:
 - 80 kbp, 70 kbp, 50 kbp, 40 kbp, 30 kbp, 20 kbp, 10kbp, 5kbp
 - Sum size of assembly B is 305 kbp
 - N50 contig length is 50 kbp

Overview

- DNA Sequencing Applications
 1. Resequencing
 2. Capture/Targeted Resequencing
 3. RAD/ddRAD/GBS
 4. Bisulfite Sequencing
 5. De novo assembly
- Reference Genome availability
- N50

Python
functions!